



Pipet Help



Powerful yet quite pump

Soft grip push buttons

Double safety valve

Autoclavable pipet cone

The smart Pipet Controller

Pipet Help

- Can fill 25 ml pipet under 4 seconds
- Safety valve and hydrophobic filters provide double protection against fluid penetration
- Two different speed modes - High & Low along with Gravity drain
- The pump speed can also be fine-tuned by varying finger pressure on operating knobs for better control of speed
- Cadmium free environment friendly NiMH batteries
- Batteries can be changed very easily by the user
- The intelligent charger prevents over charging / heating of batteries
- Low battery indicator

Technical Data

Power	2 x 1.2 V / 1000 mAh NiMH batteries
Battery Service Life	Eight hours continuous use
Battery Charging Time	2 - 3 hours to charge fully
Volumetric Range	Plastic or Glass pipette from 0.1 ml up to 200 ml
Safety System	Replaceable hydrophobic 0.45 µm PTFE filter Safety check valve
Standard Accessory	Battery charger (UK plug) Hydrophobic PTFE membrane filter 0.45 µm 1Pc Hydrophobic PTFE membrane filter 0.20 µm 1Pc Wall mounting stand with screw


Ordering Information

PipetHelp

<u>Cat#</u>	<u>Description</u>
PH01-R	PipetHelp - Red
PH01-B	PipetHelp - Blue
PH01-O	PipetHelp - Orange
PH01-G	PipetHelp - Green

Spares for PipetHelp

<u>Cat#</u>	<u>Description</u>
PH01-10	Silicone pipet adaptor with safety valve
PH01-11R	Pipet cone- autoclavable - Red
PH01-11B	Pipet cone- autoclavable - Blue
PH01-11O	Pipet cone- autoclavable - Orange
PH01-11G	Pipet cone- autoclavable - Green
PH01-12	PTFE membrane Filter Set, 0.20 µm (5 pcs)
PH01-13	PTFE membrane Filter Set, 0.45 µm (5 pcs)
PH01-14	Battery charger
PH01-15	Adaptor for charger plug
PH01-16	NiMH battery (set of two batteries)
PH01-17	Wall mounting stand



Comfortable stress-free operating knob



Wall mounting stand with fixing screw



One hand operations



Attractively color coded



www.accumaximum.com

Fine Care Biosystems

